



	Vibramax 100	Rotamax 120	Unimax 1010	Unimax 2010
Motion	vibrating	orbital	orbital	orbital
Rotation speed range	150 - 1,350 rpm	20 - 300 rpm	30 - 500 rpm	20 - 400 rpm
Operating Mode	timer or continuous	timer or continuous	timer or continuous	timer or continuous
Orbit / Stroke	3 mm	20 mm	10 mm	20 mm
Timer	Yes	-	Yes	Yes
Power input	31 W	33 W	50 W	115 W
Platform size	220 x 220 mm	220 x 220 mm	290 x 258 mm	390 x 340 mm
Max. Load	2 kg	2 kg	5 kg	10 kg
Overheat protection	self resetting	self-resetting	self-resetting	self-resetting
Protection class DIN EN 60529	IP 30	IP 30	IP 40	IP 20